

Алматы (7273)495-231	Иваново (4932)77-34-06	Магнитогорск (3519)55-03-13	Пермь (342)205-81-47	Тверь (4822)63-31-35
Ангарск (3955)60-70-56	Ижевск (3412)26-03-58	Москва (495)268-04-70	Ростов-на-Дону (863)308-18-15	Тольятти (8482)63-91-07
Архангельск (8182)63-90-72	Иркутск (395)279-98-46	Мурманск (8152)59-64-93	Рязань (4912)46-61-64	Томск (3822)98-41-53
Астрахань (8512)99-46-04	Казань (843)206-01-48	Набережные Челны (8552)20-53-41	Самара (846)206-03-16	Тула (4872)33-79-87
Барнаул (3852)73-04-60	Калининград (4012)72-03-81	Нижний Новгород (831)429-08-12	Саранск (8342)22-96-24	Тюмень (3452)66-21-18
Белгород (4722)40-23-64	Калуга (4842)92-23-67	Новокузнецк (3843)20-46-81	Санкт-Петербург (812)309-46-40	Ульяновск (8422)24-23-59
Благовещенск (4162)22-76-07	Кемерово (3842)65-04-62	Ноябрьск (3496)41-32-12	Саратов (845)249-38-78	Улан-Удэ (3012)59-97-51
Брянск (4832)59-03-52	Киров (8332)68-02-04	Новосибирск (383)227-86-73	Севастополь (8692)22-31-93	Уфа (347)229-48-12
Владивосток (423)249-28-31	Коломна (4966)23-41-49	Омск (3812)21-46-40	Симферополь (3652)67-13-56	Хабаровск (4212)92-98-04
Владикавказ (8672)28-90-48	Кострома (4942)77-07-48	Орел (4862)44-53-42	Смоленск (4812)29-41-54	Чебоксары (8352)28-53-07
Владимир (4922)49-43-18	Краснодар (861)203-40-90	Оренбург (3532)37-68-04	Сочи (862)225-72-31	Челябинск (351)202-03-61
Волгоград (844)278-03-48	Красноярск (391)204-63-61	Пенза (8412)22-31-16	Ставрополь (8652)20-65-13	Череповец (8202)49-02-64
Вологда (8172)26-41-59	Курск (4712)77-13-04	Петрозаводск (8142)55-98-37	Сургут (3462)77-98-35	Чита (3022)38-34-83
Воронеж (473)204-51-73	Курган (3522)50-90-47	Псков (8112)59-10-37	Сыктывкар (8212)25-95-17	Якутск (4112)23-90-97
Екатеринбург (343)384-55-89	Липецк (4742)52-20-81		Тамбов (4752)50-40-97	Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

www.getinge.nt-rt.ru | | gtw@nt-rt.ru

Технические характеристики на ротационные запайщики Proseal, Pack Rotary Sealer, Proseal Plus, Proseal Premium компании GETINGE

Getinge Pack Rotary Sealer

User-friendly rotary sealer supporting high-volume sealing operations

Getinge Pack Rotary Sealer can be used with sterilization packaging made of medical paper/film and Tyvek® material/film in compliance with EN 868-5 and ISO 11607-1.

Packaging small and large instruments can be time consuming and sometimes troublesome. We have designed a range of sealers that maximize both productivity and effectiveness, improving throughput at your medical facility.



Getinge Pack Rotary Sealer offers:

Safe and secure seal

- Easy sealing
- Sealed in an automated process ensuring repeatable and effective sealing results
- The pouch guide allows the user to lean the pouch for a straight and perfect alignment for entry

Technical data

Sealing Parameters	Getinge Pack Rotary Sealer 220V	Getinge Pack Rotary Sealer 110V
Temperature range (adjustable)	50°C – 250°C (122°F– 482°F)	
Speed range (adjustable)	10 m/min (33 ft/min)	
Pressure range adjustable (mechanically)	40–100 N	
Seal width	14	
Seal length	Unlimited	
Features		
Automatic seal start	Yes	
Control	Microprocessor	
Seal type	Multiline seal (4 lines)	
Heat up time to 180 °C (356 °F)	180 s	
Material	Sanded stainless steel	
Display	LCD Screen	
Power and technical information		
Electricity	200–240 V, 50–60 Hz	110–120 V, 50–60 Hz
Power	650 W	
Weight	10 kg (22 lb)	
Dimensions (W x D x H)	425 x 225 x 155 mm (16 x 9 x 6 in)	



Transfer Roller

The transfer roller is a useful accessory for Getinge Pack Rotary Sealer, providing safe and easy handling with heavy packs and during high-volume packaging operations.



Cutting Device

The roll cutting device comes in three different sizes, for wall or tabletop use. The sliding-action knife is completely built into the housing for safe handling.

Getinge Proseal

A range of rotary sealers supporting high-volume sealing operations with validatable sealing processes.



Getinge's rotary sealers provide validated sealing of sterile bags and pouches

To make it easier to seal a perfect seam Getinge has developed a slow-start function that minimizes the number of wrinkled seams. Use the slow-start function on the touchscreen and the belt slows down, making it easier to insert the pouch. Furthermore, the sealers are tilted for a better ergonomic position and have easy to read signals from LED-light and panel.

The Getinge ProSeal rotary model is available in two versions: Plus and Premium, both validated according to DIN EN ISO 11607-2

Intuitive user interface

The intuitive, user-friendly design of the 4.3" Centric interface has been developed to easily change settings.

Integrated printer

Getinge Proseal Premium models are equipped with a built-in module to print traceability information on pouches. Up to nine different parameters can be printed.

Slow-start function

A unique slow-start function automatically decreases the sealing speed. This helps to feed the heavy and bulky pouches to optimize seal quality.

T-DOC

Proseal Sealers can connect with T-DOC Getinge's traceability system. All critical parameters can be logged and monitored for traceability and documentation purposes.

All logged information can be accessed and transferred via USB or Ethernet connection.

Technical data

Sealing parameters:	Getinge Proseal Plus	Getinge Proseal Premium
Temperature range (adjustable)	100 °C–200 °C (212 °F–392 °F)	
Speed range (adjustable)	120–170 mm/s (4.7–6.7 in/s)	
Pressure range adjustable (mechanically)	40–100 N	
Seal width	15 mm (0.59 in)	
Seal length	Unlimited	
Power and technical information:		
Electricity	220–240 V, 50–60 Hz (110–127 V, 50–60 Hz)	
Power	550 W	
Weight	9.8 kg (21.6 lbs)	10.1 kg (22.3 lbs)
Dimensions (W x D x H)	580 x 255 x 133 mm (23 x 10 x 5.25 in)	580 x 255 x 155 mm (23 x 10 x 5.25 in)
Features:		
Slow-start function	Yes	
Standby time (adjustable)	1-60 min	
Heat up time to 180 °C (356 °F)	180 sec	
Automatic seal start	Yes	
Touchscreen	Capacitive scree, colored	
Control	Microprocessor	
Connectivity	USB, Ethernet	
Printer	-	Yes (Integrated)



Support table

The support table enables safe handling, even with heavy packs and high-volume packaging operations. It is especially designed for the Proseal Premium Sealer.



Cutting device

The roll cutting device comes in three different sizes, for wall or tabletop use. The sliding-action knife is completely built into the housing for safe handling.

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97

Тверь (4822)63-31-35
Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

www.getinge.nt-rt.ru | gtw@nt-rt.ru